The lot-specific document is available at the SERO website

EN

Intended use

Seronorm® Cardiac Acute Liq is intended for use as a quality control of cardiac and acute markers to monitor immunoassay measurement procedures in a laboratory environment by professional users.

Summary

Seronorm® Cardiac Acute Liq is supplied in 4 levels: L-1A, L-1B, L-2 and L-3. The product is suitable for use as a quality control for the analytes listed in the data tables in the lot-specific product documentation.

Composition

This product is a liquid human-based control serum produced from blood collected from voluntary blood donors.

Warnings and precautions

For in vitro diagnostic use by health care professionals.

Each human donor unit has been separately tested by CE-marked or FDA-approved methods to be negative for HBs antigen, HIV p24-antigen and HIV I, II and HCV antibodies. However, since no method can completely exclude the presence of infectious agents, this material should be handled as though capable of transmitting infectious diseases and disposed of accordingly. Exercise the normal precautions required for good laboratory practice.

Increased turbidity might indicate bacterial contamination, in which case the control material cannot be used and must be disposed of safely.

In case of any serious incident related to the use of the quality control, please report it immediately to SERO. Contact us at support@sero.no or call +47 66 85 89 00. Additionally, notify your competent authority.

The Safety Data Sheet (SDS) is available at our website www.sero.no/theserolibrary or on request by using the email address sales@sero.no.

Instructions for use

To secure uniformity, ensure that each vial is similarly treated during preparation for use. Always mix contents thoroughly before use.

1. Place the vial on a roller mixer for approximately 20 minutes at room temperature until the contents have thawed completely. Avoid foaming of the material during mixing. Do not shake.
2. Remove the cap and stopper from the vial.
3. Remove the amount of material required for testing. Replace the stopper and cap. Promptly store the original vial at 2-8 °C. For optimal analyte stability, analyze the material immediately.
4. Prior to subsequent use, remove the vial from 2-8 °C and mix for a few minutes before sampling. Replace the stopper and cap, and promptly store the original vial at 2-8 °C after use. For optimal analyte stability, analyze the material immediately.
5. DO NOT RETURN control material to the original vial once removed.
6. Establish target values and analytical ranges for the control material based on the measurement procedures in the lab.

Storage & Stability

Unopened vials: Store at ≤ -20 °C. Stable to expiration date as printed on the vial. Opened vials: Please refer to the lot-specific product documentation available online.

Limitations

Refer to the lot-specific product documentation for specific analyte and stability claims.

Assignment of values

Refer to the lot-specific product documentation available online for specific analyte claims.

Visit www.sero.no/theserolibrary to download the lot-specific product documentation and a copy of this document. The 8-digit number printed on the label at the top of the first page of this document is the download code required to access the online lot-specific documentation. To open the electronic documents, Adobe PDF Reader is required. This can be downloaded at no cost from get.adobe.com/reader. If you do not have internet access, the label is missing or if you are unable to download the documents from the website, please contact your SERO representative or SERO directly at sales@sero.no.

Disposal of waste

This quality control material must be disposed of in accordance with applicable national, regional or local regulations for waste. Waste should be handled with the same care as biohazardous waste material.

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